

Institute for Systems Biology Protocol
Direct Incorporation of Cy3/Cy5 During Reverse Transcription
Written by Anne Rosenwald and Todd Eckdahl; August 2003)

Materials Needed:

Isolated and quantitated total RNA samples
Microarray slides (70-mer plus-strand oligomers)
RNase-free water
oligo dT primer (16- to 18-mer) at 1 μ g/ μ l
Coverslips, 22 x 40mm size from Corning
100 mM DTT (dithiothreitol)
low dTTP dNTP mix (10 mM each dATP, dCTP, dGTP, **1 mM dTTP**)
Cy3-dUTP and Cy5-dUTP (1 mM each [separately])
3 M Ammonium Acetate, pH 5.2
100% Ethanol, 70% Ethanol

Superscript II Reverse Transcriptase, 5X first strand buffer
RNase A (4 mg/ml)
RNase H (2 unit/ μ l)

Reverse Transcription and Cy-dye Incorporation

1. In duplicate, aliquot 50 μ g total RNA (one for each treatment)
 - a. already checked with denaturing agarose gel
 - b. quantitated with UV spectrophotometer
 - c. precipitated (eg. 1/10 volume 3M NaOAc pH 5.2, 2 volumes EtOH)
2. To each tube, add 2.5 μ g oligo dT
3. Adjust volume to 11 μ l with DEPC-treated H₂O
4. Heat to 75°C for 10 min
5. Cool slowly to room temperature and spin down (Note: keep at RT from this point on)
6. Add the following in order:
 - a. 4 μ l Superscript first strand 5X buffer
 - b. 2 μ l DTT (100 mM)
 - c. 1 μ l dNTPs (10 mM each dATP, dCTP, dGTP and **1 mM dTTP**)
 - d. 1 μ l Cy-dye labeled dUTP (1 mM) (One gets Cy-3 dUTP and one gets Cy-5 dUTP)
 - e. 1 μ l Superscript Reverse Transcriptase II (200 units/ μ l; make sure this is Exonuclease-free)
7. Mix gently and incubate at room temperature for 10 min
8. Incubate at 42°C for 2-3 hours (no more than ~5 hours; do not do this step overnight)
9. Heat sample to 95°C for 2 min
10. Place samples on ice, spin down (can store at -20°C at this point if necessary)

Degrade RNA

1. Make sure contents of tubes are spun down
2. Add 0.5 μ l of RNase A (4 mg/ml) at room temperature (Promega)
3. Add 0.5 μ l of RNase H (2 U/ μ l) (Fermentas) (Note: Not clear if really necessary to use both enzymes; the RNase H is fairly expensive)
4. Incubate at 37°C for 15-30 min

(Alternative RNA Degradation)

1. Add 3.5 μ l 0.5 M NaOH/50 mM EDTA, then
2. Heat at 65 C, 10 minutes
3. Add 5 μ l 1 M Tris 7.5

Purification (Using Qiagen PCR CleanUp Kit)

1. Add 25 μ l high-quality H₂O to samples (add only 17.5 μ l if alternative RNA degradation was done)
2. Add 2.7 μ l 3 M Sodium Acetate, pH 5.2
3. Add 250 μ l QIAquick buffer PB
4. Apply each sample to a QIAquick column (the DNA should stick to the column here)
5. Centrifuge for 30 sec at full speed
6. Take the column flow-through and replace back onto the top of the column and spin a 2nd time
7. Place the flow-through back in the original tube and save in case of problems with the purification.
8. Wash with 400 μ l QIAquick buffer PE, spin 30 sec at full speed and discard flow-through (your DNA remains on the column)
9. Repeat step 8, discarding flow-through
10. Spin the column briefly once more to get rid of remainder of wash solution
11. Place column in a clean, well-labeled 1.5 ml elution tube
12. Apply 30 μ l buffer EB to center of column without touching the membrane
13. Wait one min, then centrifuge 1 min at full speed (gradually increase from 0 to full speed to avoid shearing off the eppendorf tube lids). Your DNA is in the flow-through this time!
14. Again apply 30 μ l buffer EB to center of column without touching the membrane. DNA should be in a volume of 60 μ l.

Determining the incorporation of labeled nucleotides

1. Use a 384-well spec plate (if available) so that small volumes can be used
2. Use undiluted labeled sample and buffer as the blank (~40 μ l of each). Read absorbances at 260 nm, 280 nm, 550 nm for Cy3 and 650 for Cy5.
3. Calculations:

- a. Extinction coefficients are 150,000 for Cy3 and 250,000 for Cy5
- b. Purity: Corrected A260/A280 (want this to be ~1.8)
- c. Pmol dye/ μ l of sample = A(dye)/(extinction coefficient * 10⁶)
- d. dsDNA = 50 μ g/ml
- e. ssDNA = 33 μ g/ml
- f. ssRNA = 40 μ g/ml

Example:

$$\text{Total dsDNA } (\mu\text{g}) = 50 \mu\text{g/ml} * 1 \text{ ml} / 1000 \mu\text{l} * A_{260} * \text{Volume } (\mu\text{l}) * \text{DF}$$

4. If labeling efficiencies are fairly close for each dye and the nucleic acid yield is also similar, then about 30 pmol of dye of the sample per hybridization is a good starting point to product good intensity on most arrays. It is better to control the amounts of starting RNA/cells and have a good internal control and normalization scheme to deal with incorporation discrepancies.

Institute for Systems Biology
Direct Incorporation of Cy3/Cy5 During Reverse Transcription

Cost Analysis

Reagent	Supplier	Startup Cost	Number of Two Label Expts	Cost per Two Label Expt
oligo dT	Operon	\$75 for 100 ug	20	\$3.75
Superscript II	Invitrogen Life Technologies	\$220 for 10,000 units	25	\$8.80
Cy3-dUTP	Amersham	25 nmoles for \$295	25	\$11.80
Cy5-dUTP	Biosciences	25 nmoles for \$304	25	\$12.16
Rnase A	Promega	\$212 - 1 ml, 4 mg/ml	1000	\$0.22
Rnase H	Promega	\$146 for 50 units	25	\$5.84
QIAquick	Qiagen	50 columns for \$76	25	\$3.04
Totals		\$1328 for 25 Experiments		\$45.61 per Expt